

GD05

Natural latex powder free disposable glove



Sectors:

- Automotive •Construction & Utility •Engineering and Manufacturing •Food Manufacturing and Preparation •Janitorial •Hair and Beauty

Features

Manufactured from: Natural rubber latex

Flexible material: For improved feel and fit reducing fatigue during extended wear

Chlorinated: To reduce latex protein content

Food Contact: Tested in accordance with European standard EN1186 making it suitable for food contact

Powder free: Reduces the risk of dust contamination

Textured finish: For improved grip

Pat Coakley Medical 



For further technical information, please contact Customer Services on 033 33208550 or the Technical team on 033 33208555

Polyco Healthline, South Fen Business Park, South Fen Road, Bourne, Lincolnshire, PE10 0DN

GD05

Code	GD05
Colour	Natural
Sizes	XS-XL



Inners dimensions (mm)
230 x 120 x 65mm

Outers dimensions (mm)
335 x 250 x 238mm

Size	Inner Barcode	Outer Barcode
XXS/6	5024951004370	5024951005629
S/7	5024951003960	5024951003977
M/8	5024951003984	5024951003991
L/9	5024951004004	5024951004011
XL/10	5024951004028	5024951004028

Standards

Contains accelerators and latex which may cause allergic reactions, including anaphylactic response

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Product Food Approved to EC1935/2004

EU Type Examination carried out by SATRA

Tested to European Standards EN420:2003+A1:2009 & EN ISO 374-1:2016 & EN ISO 374-5:2016

Pat Coakley Medical Ltd.

Units E3 E4 Network Enterprise Park
Kilcoole A63 H244 Co Wicklow
IRELAND

+353 12018070
sales@patcoakleymedical.ie

Packing

100 gloves per dispenser
10 dispensers per case
84 cases per pallet

