



L'ORIGINALE DAL 1983
THE ORIGINAL SINCE 1983

COLLECTION





REBEL

A NEW VERSATILE AND ADJUSTABLE CLOG

Calzuro *REBEL*[®], an attractive look for a new concept of professional clogs:
the safety and comfort of Calzuro, the aerodynamic design of a sports footwear.



CE CERTIFIED



ECO-FRIENDLY



100% RECYCLABLE



LATEX-FREE



SLIP-RESISTANT



AUTOCLAVABLE



ANTI-STATIC



CRUELTY-FREE



MADE IN ITALY



MACHINE WASH UP
TO 100°C / 212°F



UNISEX



- AUTOCLAVABLE UP TO 134°C / 273°F
- LATEX-FREE
- LIGHT
- 100% RECYCLABLE
- MADE IN ITALY
- MEDICAL DEVICE CLASS I
- NON STERILE: DM 2431858



SLIP-RESISTANT SOLE

According to
EN ISO 20347:2012
SRC requirements

ANTISTATIC
PLUG
on the sole



Use the heel tab
positioned on
the T-CLOSE®
to put on and
take off the clogs

T-CLOSE®

The patented back closure for the protection of the heel, adjustable and not removable, makes the clog more stable and safer. Strap with 3 levels of regulation. Provided with heel tab to best put on and take off the clogs.



COMFORT INSOLE, REMOVAL AND ANTISTATIC


- Anatomical shape, with CUSHION ENERGY insert in heel region for cushion effect
- Arch support to maintain the correct foot position
- Designed to reduce plantar and muscular fatigue, and metatarsalgia disorders
- The upper fabric and the thermal comfort foam help maintain dry and cool foot
- Hand-wash with neutral soaps
- Entirely antistatic



FOR A
BETTER
POSTURE

AUTOCLAVABLE INSOLE, REMOVABLE AND ANTISTATIC

- Anatomical shape for greater comfort
- Autoclavable up to 134°C / 273°F (insole must remain inside the clog during autoclave cycle)
- Provided with massaging bumps to help microcirculation
- Permeable thanks to the small holes in the front
- Shock absorbing
- With antistatic plug





AVAILABLE COLORS



WHITE / BLUE



BLACK / GREEN



PINK / PURPLE
available up to 40/41



METALLIC BLUE / WHITE



GREEN / PASTEL GREEN

AVAILABLE COMBINED SIZES

36/37 – 37/38 – 38/39 – 39/40 – 40/41 – 42/43 – 44/45

CALZURO® Classic



Professional clogs made of a thermo-techno polymer compound latex-free, completely antistatic, with side holes for foot perspiration, available with or without upper holes, autoclavable to 134°C / 273°F, washable to 90°C / 194°F. Slip-resistant sole, CE certified. All colors are available with heel strap, assembled on demand. 100% Made in Italy.

Certified according to EN ISO 20347:2012 requirements:

- **A** (antistatic properties)
- **E** (energy absorption of seat region)
- **OB** (basic requirements)
- **SRC** (slip resistance)

MEDICAL DEVICE CLASS I NON STERILE: DM 2431830

SIZES

34/35 – 35/36 – 36/37
 37/38 – 38/39 – 39/40
 40/41 – 42/43 – 44/45
 46/47

**ALL COLORS ARE AVAILABLE
 WITH HEEL STRAP ASSEMBLED
 ON DEMAND**

COLORS

COLORS AVAILABLE
 WITH OR WITHOUT UPPER HOLES



WHITE



GREEN



LIGHT BLUE



BLACK



PINK

AVAILABLE UP TO 42/43



DARK BLUE



MET. BLUE



RED



MAROON



MET. GREY

COLORS AVAILABLE
 WITHOUT UPPER HOLES ONLY



PASTEL
 PINK
 AVAILABLE UP
 TO 40/41



PASTEL
 GREEN
 AVAILABLE UP
 TO 42/43



YELLOW



ORANGE



PURPLE

COLORS AVAILABLE
 WITH UPPER HOLES ONLY



TEAL



FANCY
 WHITE-BLUE



FANCY
 WHITE-PINK
 AVAILABLE UP TO 42/43



CE CERTIFIED



ECO-FRIENDLY



100% RECYCLABLE



LATEX-FREE



ANTISLIP



AUTOC LAVABLE



ANTISTATIC



CRUELTY-FREE



MADE IN ITALY



MACHINE WASH
 UP TO 90°C/194°F



UNISEX

ACCESSORIES

Classic

INSOLE



Arch support insole for Calzuro Classic in special foam that allows feet to breathe, autoclavable to 121°C / 250°F up to 50 times.

SIZES

34 - 35 - 36 - 37 - 38 - 39 - 40 - 42 - 44 - 46




CALZURELLI




SWAROVSKI

HEEL STRAP



COLORS



NB: FANCY COLORS SUPPLIED ONLY WITH WHITE HEEL STRAP
AUTOCLAVABLE TO 134°C / 273°F

CALZURO®
Fancy



Professional clogs made of a thermo-techno polymer compound latex-free, completely antistatic, with side holes for foot perspiration, available with upper holes only, autoclavable to 134°C / 273°F, washable to 90°C / 194°F. Slip-resistant sole, CE certified. Available with heel strap, assembled on demand. 100% Made in Italy.

Certified according to EN ISO 20347:2012, requirements:

- **A** (antistatic properties)
- **E** (energy absorption of seat region)
- **OB** (basic requirements)
- **SRC** (slip resistance)

MEDICAL DEVICE CLASS I NON STERILE: DM 2431830

SIZES

34/35 – 35/36 – 36/37
37/38 – 38/39 – 39/40
40/41 – 42/43 – 44/45
46/47

COLORS



FANCY
WHITE-BLUE



FANCY
WHITE-PINK
AVAILABLE UP TO 42/43

AVAILABLE WITH WHITE HEEL STRAP ASSEMBLED ON DEMAND



CE CERTIFIED



ECO-FRIENDLY



100% RECYCLABLE



LATEX-FREE



ANTISLIP



AUTOC LAVABLE



ANTISTATIC



CRUELTY-FREE



MADE IN ITALY



MACHINE WASH
UP TO 90°C/194°F



UNISEX

CALZURO® Light



Professional clogs made of ultra-light soft plastic latex-free material; casein-resistant, with antistatic plug and side holes for foot ventilation. Available only without upper holes and with heel strap. Machine washable up to 50°C / 122°F. Slip-resistant sole, CE certified. 100% Made in Italy.

Certified according to EN ISO 20347:2012, requirements:

- **A** (antistatic properties)
- **E** (energy absorption of seat region)
- **OB** (basic requirements)
- **SRC** (slip resistance)

MEDICAL DEVICE CLASS I NON STERILE: DM 2431857

SIZES

35 - 36 - 37 - 38

39 - 40 - 41 - 42

43 - 44 - 45 - 46

COLORS



WHITE



GREEN



NAVY BLUE



BLACK



PINK

AVAILABLE
UP TO 41



LILAC

AVAILABLE
UP TO 41



LIGHT BLUE



CE CERTIFIED



ECO-FRIENDLY



LATEX-FREE



ANTISLIP



CASEIN-RESISTANT



ANTISTATIC



CRUELTY-FREE



MADE IN ITALY



MACHINE WASH
UP TO 50°C/122°F



UNISEX

ACCESSORIES

Light



ANTISTATIC INSOLE

Comfortable antiperspirant and antistatic insole, made in recycled latex.

To be trimmed for custom feet.

Drytech fabric that helps keeping the right temperature of the foot without blocking its natural perspiration. Washable up to 30°C / 86°F with mild soap.

CALZURO®

Aqua

Flip flops made of EVARITE®, a new, soft, lightweight, latex-free plastic material which does not favor fungi and bacteria proliferation. Anatomically shaped footbed, well-dressing anatomic shape, slightly raised heel and toe bed. Machine washable up to 50°C / 122°F. Slip-resistant sole, CE certified. 100% Made in Italy.

Certified according to EN ISO 20347:2012, requirements:

- **E** (energy absorption of seat region)
- **OB** (basic requirements)
- **SRC** (slip resistance)

MEDICAL DEVICE CLASS I NON STERILE: DM 1831009

SIZES

35 – 36 – 37 – 38 – 39 – 40
41 – 42 – 43 – 44 – 45 – 46

COLORS



CE CERTIFIED



ECO-FRIENDLY



LATEX-FREE



ANTISLIP



CRUELTY-FREE



MADE IN ITALY



MACHINE WASH
UP TO 50°C/122°F



UNISEX

CALZURO® **Pro Safety**

WITH SAFETY TOE-CAP



Soft and practical professional footwear with closed back-heel, made of a brand new special ultra-light copolymer, latex-free and casein-resistant.

Easy wear slip on and slip off and safe fit for maximum comfort and further your safety policy at work. Supplied without side vents, with antistatic plug and with synthetic polymer metal-free protection toe-cap (impact resistance of the toe-cap: 200 Joule on the toe).

The safety toe-cap provides adequate protection against mechanical risks such as impact and crashing/breakthrough, at the same time it is lighter than leather safety shoes with steel toe-cap.

Slip-resistant sole, CE certified. 100% Made in Italy.

Certified according to EN ISO 20345:2011, requirements:

S4 (additional requirements):

- **Closed heel area**
- **A** (antistatic properties)
- **E** (energy absorption of seat region)
- **Cleated outsole**
- **SRC** (slip resistance)

MEDICAL DEVICE CLASS I NON STERILE: DM 1770386

SIZES

36 – 37 – 38 – 39 – 40

41 – 42 – 43 – 44 – 45

46

COLORS



BLACK



CE CERTIFIED



ECO-FRIENDLY



LATEX-FREE



ANTISLIP



CASEIN-RESISTANT



METAL-FREE
SAFETY TOE-CAP



ANTISTATIC



CRUELTY-FREE



MADE IN ITALY



MACHINE WASH
UP TO 50°C/122°F



UNISEX



EASY-FIT

CALZURO® Effetto

Slippers for medical and paramedical personnel, slip-resistant PU (polyurethane) sole, leather upper made of pigmented perforated leather for a healthy and natural foot perspiration, soft padding on the instep. Can be supplied with heel strap.

Certified according to EN ISO 20347:2012, requirements:

- **A** (antistatic properties)
- **E** (energy absorption of seat region)
- **SRC** (slip resistance)

MODELS



WITHOUT HEEL STRAP



WITH HEEL STRAP



SIZES

35 – 36 – 37 – 38 – 39 – 40
41 – 42 – 43 – 44 – 45 – 46

COLORS



WHITE
(with white sole)



BLUE
(with black sole)



ANTISTATIC



UNISEX



MADE IN ITALY



CE CERTIFIED



ANTISLIP

CALZURO®

Bab

Antistatic white boots in special polymer compound (nitrile rubber), non-toxic, resistant to hydrocarbons (fuel oil). A good washing resistance at medium-high temperatures 85°C / 185°F helps maintain within optimal limits the bacteria spread. Resistant to animal fats, vegetable oils and casein, slip-resistant sole, knee-high. White polyester knitting inside.

Washable up to 85°C / 185°F. 100% Made in Italy.

Certified according to EN ISO 20347:2012, requirements:

- **O4** (closed seat region, antistatic properties, energy absorption of seat region)
- **FO** (bottom sole resistance to fuel oil)
- **SRC** (slip resistance)

SIZES

36 – 37 – 38 – 39 – 40 – 41 – 42

43 – 44 – 45 – 46 – 47 – 48

COLORS



WHITE



CE CERTIFIED



ANTISLIP



CASEIN-RESISTANT



ANTISTATIC



CRUELTY-FREE



UNISEX



MACHINE WASH
UP TO 85°C/185°F



MADE IN ITALY



FUEL OIL - RESISTANT



HOW TO SANITIZE CALZURO® CLOGS





REBEL

- Sterilizable in autoclave up to 134°C / 273°F (except comfort removable insole)
- Washable in washing machine up to 90°C / 194°F (except comfort removable insole)
- Washable in dish washer up to 70°C / 158°F (except comfort removable insole)



CLASSIC MODELS

CLASSIC with or without upper holes, FANCY, PASTEL

- Sterilizable in autoclave up to 134°C / 273°F
- Washable in clogwasher up to 90°C / 194°F with specific disinfecting cleaner
- For immersion in 1% sodium hypochlorite aqueous solution for 15 minutes and subsequent rinsing
- By spraying of pure alcohol at 60-80% and leaving it for at least 15 minutes (optional subsequent rinsing)
- For Immersion with specific disinfectant (such as: Napisan®), thoroughly following the manufacturer's instruction and timing
- Hand washing at 30°C / 86°F minimum for at least 3 minutes with whichever soap and leaving it in for at least 15 minutes before rinsing



MODELS IN EXPANDED EVA

LIGHT – PRO SAFETY – AQUA

- Washable in clogwasher up to 50°C / 122°F with specific disinfecting cleaner
- For immersion in 1% sodium hypochlorite aqueous solution for 15 minutes and subsequent rinsing
- By spraying of pure alcohol at 60-80% and leaving it for at least 15 minutes (optional subsequent rinsing)
- For Immersion with specific disinfectant (such as: Napisan®), thoroughly following the manufacturer's instruction and timing
- Hand washing at 30°C minimum for at least 3 minutes with whichever soap and leaving it in for at least 15 minutes before rinsing





OUR PRODUCTION, OUR PRIDE

Made in Italy <

Updated CE certifications <

High Quality Materials <

Designed for every need <

For over 40 years we have devoted ourselves to every phase of the manufacturing - entirely MADE IN ITALY - with great attention. From the selection of new materials to the design and choice of avant-garde colors, we carefully follow every step of the way. We are constantly up to date on European regulations and certifications, in order to offer high quality and maximum safety.

The passion, attention to details and innovation lead our work; and we proudly value the excellence we have developed in these 40 years. This expertise allows us to be at the side of professionals step by step and to offer high quality solutions for every need.

These are the reasons why Calzuro has always been the reference brand for professional footwear.

**FET
GOOD**





L'ORIGINALE DAL 1983
THE ORIGINAL SINCE 1983

100% MADE IN ITALY
www.calzuro.it

BIHOS S.R.L.
Via Dell'Artigianato, 18
35010 VIGONZA (PD)

C.F. e P.IVA 01865630287
Tel. 049-8930259

www.bihos.com
bihos@bihos.com
info@bihos.com

